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The Beconfident® Multifunctional formula ensures clean, visibly whiter teeth, strengthens the enamel, freshens the breath and protects against plaque and cavities.

The various flavours + Mint flavour relaxes and recharges mind and body.

How to use: Brush the teeth daily in the morning and in the evening for at least 2 minutes with the Beconfident® Multifunctional Whitening Toothpaste. Add a sip of water and swirl around for additionally 30 seconds. Best feeling ever!

b beconfident®

Beconfident® – Oral care innovation

Group Beconfident Europe AB

Tennisvägen 31, SE-269 33 BÅSTAD, Sweden

www.beconfident.com

b beconfident

Sonic Silver

Whitening Toothbrush

Quick Guide

Congratulations on your quality purchase from Beconfident® – We hope you will love our Sonic toothbrush and toothpastes

Before using this product, read this user manual carefully to learn about the features and the recommendations (the do's and don'ts) and how you can make the best use of it. Note: Charge the Beconfident sonic toothbrush completely prior to using it for the first time.

EN Quick guide how to use the brush:

1. Press power button to turn on the brush
2. Press again in less than 5 seconds to toggle through brushing modes.
3. Long press power button for 2 seconds to turn off the brush after your favourite mode is selected.
4. Wet the bristles and apply a small amount of toothpaste.
5. Place the toothbrush bristles against the teeth and (do not press hard).

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