

Teeth Whitening Mouth Trays

Art no: 1213xx & 1214xx

Individually moldable mouth trays for perfect fit

Claims:

- Works with any type of gel
- Shapes in hot water for individual fit
- Can be reused 100+ times
- With gel applicator

Product description:

Beconfident's original thermoform mouth trays are simple and easy to use. Contains two mouth trays and one gel applicator. Place them in warm water and then mold them to the shape of your teeth, giving you an accurate fit when completing the whitening process. The mouth guard can also be used as a dental guard to prevent teeth from grinding at night.



Working mechanism in short:

For the gel to be the most effective, it is important it has a close contact to the teeth surface. This is why it's important to shape the trays using hot water (approximately 80 degrees) for an individual and close fit. Once the trays are shaped, they can be stored and reused 100+ times, with any available gel.

Material:

The trays are made of soft EVA (Ethylene-Vinyl Copolymer). This product is clinically and dermatologically tested and approved not be harmful for public health or the environment. This is a CE approved product and made in a FDA approved factory that follows GMP of ISO22716:2007 Cosmetic guidelines.

The Mouth Tray Kit contains the following items:

- 1) 2 x Beconfident Original Mouthtrays
- 2) 1 x Gel Applicator
- 3) 1 x multilingual manual
- 4) 1 x link to video instructions

Beconfident is not conducting any testing on animals

Product Data:

Art No	Dimension HxWxD (mm)	Gross Weight	Net Weight	Units Display Box	Units mastercart.	Unit EAN
1213xx & 121497	175x120x33	44g	15g	8	128	7350064167045
CPNP Code	GMP ISO 22716	Country of Origin	Quality Test Pharmacies	Box bio-degradable	Pallet	RRP Euro incl VAT(+/-)
	Follows			Yes	12	15,90 €uro

Group Beconfident Europe AB

Tennisvägen 31, SE-26933 Båstad, Sweden

E-mail: info@beconfident.info Phone: +46-(0)702 56 50 00 VAT no: SE 556853-0041-01

Corporate site: www.beconfident.info, www.beconfident.co.uk, www.beconfident.se & many more

Beconfident- Top Quality Teeth whitening and Skincare systems are designed in Sweden & Made in USA and Sweden.

Beconfident® + WhiteAmin® + SimpleSmile® are registered trademarks in Europe and Beconfident® is a registered trademark in United States of America. IOS and Android Intelligent Apps belong to Beconfident Skincare AB. Other brands are pending.



Group Beconfident Europe AB

Tennisvägen 31, SE-26933 Båstad, Sweden

E-mail: info@beconfident.info Phone: +46-(0)702 56 50 00 VAT no: SE 556853-0041-01

Corporate site: www.beconfident.info, www.beconfident.co.uk, www.beconfident.se & many more

Beconfident- Top Quality Teeth whitening and Skincare systems are designed in Sweden & Made in USA and Sweden.
Beconfident® + WhiteAmin® + SimpleSmile® are registered trademarks in Europe and Beconfident® is a registered trademark in United States of America. IOS and Android Intelligent Apps belong to Beconfident Skincare AB. Other brands are pending.